



# ROMA™ MATTRESS COVER – STANDARD / COTTON EDGE / NIMBUS

This guide covers all products with the prefix RRMCM in the product code.

With or without assistance this product makes moving and/or repositioning in bed that much easier. The user can lie on the **Roma™ Mattress Cover (Roma™)** directly or use with a bed sheet placed over it. Model shown below is **Roma™** with Cotton Edge on a divan bed but they fit standard hospital beds also. The techniques shown apply to all types of our **Roma™**.

## USER INSTRUCTIONS



# 1

### PERSON UNASSISTED

Remove the **Roma™** from its packaging and place over the bed mattress, as shown.

The user can now sit on the edge of the bed, prior to rotating on the **Roma™**.



The person can now reposition themselves accordingly and sleep comfortably without risk of tissue damage or sweating. The **Roma™** fabric enables them to turn more easily.

# 2



# 3

### PERSON ASSISTED (a Flat Slide Sheet is also required)

With the **Roma™** placed over the mattress, carers can reposition the person by inserting a **Flat Slide Sheet**, as shown (on top of the **Roma™** but under the patient).

If lying on a bed sheet, insert the **Flat Slide Sheet** between the **Roma™** and bed sheet and remove after repositioning.

Unroll the **Flat Slide Sheet** under the patient without disturbing them. The **Flat Slide Sheet** can then be used to reposition the person sliding on the **Roma™** below.

# 4



# 5

Use the **Flat Slide Sheet** to reposition easily.

Then, when the person is in the desired position the **Flat Slide Sheet** removes effortlessly...



...leaving the person in exactly the position you want them to be in. Alternatively use under a standard bed sheet and use a **Flat Slide Sheet** under the bed sheet for all repositioning needs. Remove the **Flat Slide Sheet** once the person is in the desired position.

# 5



**ROMA™ MATTRESS COVER STANDARD**

Fitted bed cover made entirely from Roma™ fabric.




**ROMA™ MATTRESS COVER COTTON EDGE**

The Roma™ fabric is central to the product with fitted cotton edges for additional safety at the edge of the bed.



**ROMA™ MATTRESS COVER NIMBUS**

Designed to fit over an air flow mattress. The elastic straps secure the Roma™ to the bed and does not constrict the air system at all.



**MATERIAL CONTENT**

**Standard:** 100% Polyester. **Cotton Edge:** 100% Polyester Side Panels: Polycotton. **Nimbus:** 100% Polyester with strong elastic straps.

**STORAGE AND HANDLING**

- Store in a clean dry environment that is not subject to excessive heat.
- Keep away from sunlight.
- Store flat and avoid placing heavy objects on top of the product.

**WARNINGS AND PRECAUTIONS**

- Do not leave or drop on the floor; it is a slip hazard.
- Test for slippiness before use. If over-laundered and doesn't slide please discard or replace.
- Inspect product before use for signs of damage or tears.

**CLEANING INSTRUCTIONS**

- Wash as per laundry instructions.
- Can be wiped clean, between use with the same patient, with natural water and detergent or wipes that don't contain alcohol, solvents, bleaching or abrasive agents.
- Cleaning materials used should be patient safe and biodegradable.

**DISPOSAL INSTRUCTIONS**

- Non-recyclable.
- Dispose of as clinical waste if it has been used by an infectious person.
- Otherwise dispose of through normal waste management.



Keep dry




Keep away from sunlight



Caution



Do not use if package is damaged



Storage Temperature Limit



Wash at 30-74°C. Tumble dry on low heat. Do not iron. Do not use bleach. Do not use fabric conditioner.



Patient Information website



Medical Device



**MANUFACTURER**

Hospital Direct (Marketing) Ltd, Units 2, 3 & 4, The Green, Clun, Shropshire SY7 8LG, United Kingdom

**EU REPRESENTATIVE:** Obelis S.A., Bd Général Wahis 53, 1030 Brussels, Belgium



[www.hospitaldirect.co.uk](http://www.hospitaldirect.co.uk)

**HOSPITAL DIRECT (HD) ADDITIONAL SYMBOLS**

At HD we understand that not everyone is familiar with all equipment, so to help therapists and healthcare professionals in assessing a product's suitability for a client, we have included two assessment tools, Functional Independence Measure (FIM) and Easy Guide Symbols. These will enable the healthcare professional to assess their client and decide whether or not a product will help the person based on their level of dependence and functional ability.

**FUNCTIONAL INDEPENDENCE MEASURE (FIM) ASSESSMENT**

Our FIM assessment guide for each product will help you decide the suitability of this product for the person's ability and need. Based on the standard criteria from Level 6 where the person can use the product unaided and unsupervised to Level 1 where all the assistance is provided by the carer and the client can do nothing, this guide is easy, quick and simple to use to check suitability against ability and circumstance.

**FIM LEVELS QUICK GUIDE. SUBJECT VS CARER HELP**



**EASY GUIDE SYMBOLS**

We use four symbols to indicate where a product is suitable and safe for the person to use unassisted, with a single carer or with multiple carers. These symbols indicate the minimum recommendation.



HD has a policy of continuous development and, as such, reserves the right to alter specifications (including measurements, materials and colours) without prior notice. If any serious incident has occurred in relation to the device it should be reported to the manufacturer and the competent authority in which the user and/or patient is established.